

Autosampler Vials

C0000008 C0000013

Features

- Wide opening 9-425 vials can be used on all common autosamplers, perfectly compatible with Agilent, Thermo Scientific, Waters, Varian, etc.
- Wide opening 12x32mm, Polypropylene vials have chemical resistance, which are used for atomic absorption, water and protein analysis, capillary electrophoresis, ion chromatography applications etc.



C0000013



Specifications

Cat.No.	C0000008	C0000013
Volume	2ml	2ml
Vial Size	12*32mm	12*32mm
Material	Glass	Glass
Neck Type	Screw thread	Screw thread
Vial Color	Clear	Amber
Bottom	Flat	Flat
Application	For LC/GC	For LC/GC
Package	100pcs/pk.	100pcs/pk.

Autosampler Closures

C0000143 C0000148

Features

- The closures have a fine structure and the threads are evenly distributed for easy operation.
- Pre-slit designs are available for PTFE/silicone slit shapes.



C0000143



Specifications

Cat.No.	C0000143	C0000148
Type	Open top ribbed screw cap	Open top ribbed screw cap
Cap Color	Blue	Blue
Cap Thread	9-425	9-425
Septa	Red PTFE/White silicone septa, 1mm thickness	Red PTFE/White silicone septa, 1mm thickness
Slit Shape	/	Pre-slit
Application	For LC/GC	For LC/GC
Package	100pcs/pk.	100pcs/pk.