

# Tissue Slide Dryer

TSW-1

## Features

1

The baking table is equipped with an independent high-precision temperature controller, which can directly detect the surface temperature of the baking table to ensure accurate and stable temperature control.

2

The baking table is equipped with a wax sealing table, which supports timely and convenient wax sealing protection after the wax block is sliced.

3

High throughput and high thermal conductivity separate baking table, can carry 100 slides at the same time.

4

Intelligent PID controlled heating mode, equipped with high-precision temperature sensor, and triple temperature control protection device.

5

Touch button interface, LED digital display, equipped with heating status indicator light.

6

The flexible independent split design allows the products to be freely connected and combined, saving power ports and space.



## Specification

Model	TSW-1
Baking Table Dimension	256*188*52mm
Baking Table Capacity	100 slides (10 rows, 10 slides per row)
Wax Sealing Table Dimension	74*46*5mm
Temperature Setting Range	0~99 °C
Slope Angle of Baking Table	60°
Electricity	AC 110V~120V/AC220V~240V, 50Hz/60Hz
Max. Power	350W
External Dimension	330*375*130mm
Net Weight	6kg

# Tissue Slide Dryer (Oven)

TSW-3

## Features

1

Imported single-chip microprocessor control, easy operation and reliable performance.

4

Imported pure green digital module display with clear display.

5

It is made of special black material, which is resistant to corrosion and wear.

6

The large serrated baking table is easy to put and take out, and the baking effect is excellent.

7

The new type of heating element has fast heating, uniform heating and long service life.

2

American DALLAS company temperature measurement integrated block, PID control heating mode, high precision.

3

All parameters of setting temperature, working temperature and working state are displayed in real time, and the data is more intuitive.



## Specification

Model	TSW-3
Temperature range	RT~90 °C
Temperature control precision	±1°C
Table size	340×260 mm (50 slides)
Electricity	AC 100V~120V/AC 220V~240V
Power frequency	50Hz/60Hz
Power	200W
Dimensions((W×D×H))	350×340×130mm
Weight	6kg

# Tissue Slide dryer

THP-1135



## Features



Heat conduct fast, durable.



Self memory of temperature.



Display presetting temperature & current temperature.



Corrosion-resistance of the bath.

## Description

- 1 Adopt new type of heating body, which has the advantages of fast heating conduct, long life, energy saving, etc.
- 2 With self memory of temperature function: The set temperature is automatically retained after running
- 3 Display presetting temperature & current temperature.
- 4 Made of special materials, corrosion resistance and wear resistance

## Specifications

Model	THP-1135
Temp. setting	30~90°C (constant adjusting)
Temp. accuracy	±2°C
Ambient temperature	0°C~40°C
Electricity	AC220V 50Hz, 250W
Dimension (W*D*H)	450 x 400 x 200 mm
Weight	7 kg