

Infitek

UV TRANSILLUMINATOR



Infitek



Infitek



Infitek



Infitek



Infitek

Infitek

Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: infitek.com

Email: info@infitek.com

Service: support@infitek.com

Address: Room 201, Building A, No.1 Qianwan 1st Road, Qianhai
Shenzhen-Hong Kong Cooperation Zone, Shenzhen (Settled in
Shenzhen Qianhai Business Secretary Co., Ltd.)

US Office

INFITEK INC.

522W RIVERSIDE AVE STE N,SPOKANE,WA 99201

EMAIL: INFO@INFITEK.COM



UV Transilluminator

GEP-UV1 GEP-UV2



Application

- It is suitable for excitation and recovery of nucleic acid and protein gel dyes such as EtBr, SYBR® Safe, SYBR Gold, SYBR® Green I & II, SYPRO® Ruby, SYPRO® Orange, Coomassie Fluor™ Orange stains, Gel Green, Gel Red and Lumitein™.

Features



Strong ultraviolet light excites the sample, no darkroom is required, the sample is clearly visible, and the operation is simple and convenient.



The damping 180° flat-opening ultra-strong UV protective cover completely isolates the damage of ultraviolet rays to the human.



The pressure type waterproof design prevents the buffer liquid from corroding and leaking.

Specification

Model	GEP-UV1	GEP-UV2
UV Wavelength	302nm	302+365nm
Exposure Area	21*21cm	
UV Protection	180°	
Humidity	≤90%	
Electricity	110V-220V,50-60Hz	
Dimension (L*W*H)	355*315*108mm	
N.W./G.W.(kg)	3.9kg/7.25kg	3.9kg/8kg
Package Dimension(W*D*H)(mm)	430*440*240mm	440*440*240mm