

Blood Collection Monitor

BCM-12A    BCM-12B



BCM-12A



BCM-12B

Feature



Preset volume of blood collection function  
Raring and calibration function,  
ensure measuring accuracy



If collection volume is preset, the balance  
will stop shaking when there is 10ml gap  
and get into accurate measuring status



Shaking auto stops when bags are taken  
away during collection and restart working  
after bags are remounted



Automatic alarm when reaching  
the preset value



Description

Smart Weighing Instrument for Blood Collection has the functions of load bearing , detection of bag weight , volume control tray vibration and auto restart, alarm and auto approaching- volume stop . It is suitable for all the blood collection organizations and is a convenient and accurate instrument for blood collection with plastic bags.



Unique modern sensor  
and LED display



Soft shaking mode,  
improving the survival of  
effective components in  
the blood



Humanization design  
and reasonable  
structure is convenient  
for operation

Technical Parameters

Model	BCM-12A	BCM-12B
Measuring Range	0~1200ml	
Measuring Accuracy	2%	
Display Accuracy	1ml	
Liquid Proportion	1.05g/ml	
Collection Speed	0.5~3ml/s	
Swing Frequency	33±2	
Swing Angle	13 ±2°	
Alarm	Audio and visual alarm	Audio and visual alarm with interrupter
Interrupter Thrust	/	≤14N
Power	12W	40W
Electricity	AC 110/220V ±10%, 50/60Hz	
External Size(W*D*H)	218*172*160mm	275*230*210mm
Package Dimension(W*D*H)(mm)	350*300*260mm	370*330*300mm
Net/G.W.(kg)	2.6kg/4.5kg	4.2kg/6kg