

Infitek

ANAEROBIC JAR



Infitek



Infitek



Infitek



Infitek



Infitek

Infitek

Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: infitek.com

Email: info@infitek.com

Service: support@infitek.com

Address: Room 201, Building A, No.1 Qianwan 1st Road, Qianhai
Shenzhen-Hong Kong Cooperation Zone, Shenzhen (Settled in
Shenzhen Qianhai Business Secretary Co., Ltd.)

US Office

INFITEK INC.

522W RIVERSIDE AVE STE N,SPOKANE,WA 99201

EMAIL: INFO@INFITEK.COM



PVC Anaerobic Jar, ANJ Series

ANJ125-210-3



Technical features of PVC anaerobic culture tank:



Tank body: made of transparent high-temperature-resistant polypropylene material by one-time injection molding, with high strength, high temperature resistance and strong heat permeability;



Temperature resistance: it can withstand high temperature of no more than 100 ° C
Sealing: good sealing performance, anaerobic holding time not less than 120 hours;



lock catch: Lock: unique four-button seal design, simple and convenient to operate, lock catch can be replaced separately, low maintenance cost;



Tank lid: equipped with 3 standard joint holes. It can be equipped with joints for connecting various anaerobic culture systems, and can also be used as a vent hole.



Specifications

Product Number	Capacity	Material	Applicable Specification	Standard Configuration	Optional accessories
ANJ125-210-3	Φ125/h210/2.5L	pvc	Place 1 stack of 12 dishes(Φ9-10cm) or One 300ml conical bottle	Tank body + 1 stainless steel bracket	AJ9023、A05077、AY6907、AN10

Optional accessories

Model	Name	Applicable Specification	Picture
A05077	Connector	1 Suitable for the DWSA05050 connector	
AJ9023	Connector	Suitable for the MARKI connector	
AY6907	Pressure gauge	-	
AN10	Catalyst	10g	
AN20	Catalyst	20g	
AN30	Catalyst	30g	
AN400	Catalyst	400g	

PMMA Anaerobic Jar, ANJ Series

ANJ180-240-6
ANJ130-240-3

ANJ230-240-10
ANJ130-290-3

ANJ130-120-2



ANJ130-120-2



ANJ230-240-10

Technical features of PMMA anaerobic culture tank:

- ◆ Material: The tank body and lid are made of transparent high-strength polyester material, with high smoothness, earthquake resistance and impact resistance
- ◆ Temperature resistance: It can withstand temperature of 50° C.
- ◆ Sealing: fixed and locked by fixture to fully ensure sealing.

Specifications

Model	Capacity	Material	Applicable Specification	Standard Configuration	Optional accessories
ANJ180-240-6	Φ180/h240/6L	PMMA	Place 1 stack of 12 dishes(Φ9-10cm) or four 300ml conical bottles	Tank body + 1 stainless steel bracket	AJ9023、A05077、AY6907、AN20
ANJ230-240-10	Φ230/h240/10L	PMMA	Place 3 stacks of 12 dishes(Φ9-10cm) or Eight 250ml conical bottles	Tank body + 3 stainless steel brackets	AJ9023、A05077、AY6907、AN30
ANJ130-120-2	Φ130/h120/1.35L	PMMA	Place 1 stack of 6 dishes(Φ9-10cm)	Tank body + 1 stainless steel bracket	AJ9023、A05077、AY6907、AN10
ANJ130-240-3	Φ130/h240/2.75L	PMMA	Place 1 stack of 12 dishes(Φ9-10cm)	Tank body + 1 stainless steel bracket	AJ9023、A05077、AY6907、AN10
ANJ130-290-3	Φ130/h290/3L	PMMA	Four enzyme labeled plates or cell culture plates can be placed Six identification strip culture plates	Tank body + 1 stainless steel bracket	AJ9023、A05077、AY6907、AN10



ANJ130-240-3
ANJ130-290-3



ANJ180-240-6

Optional accessories

Model	Name	Applicable Specification	Picture
A05077	Connector	1 Suitable for the DWSA05050 connector	
AJ9023	Connector	Suitable for the MARKI connector	
AY6907	Pressure gauge	-	
AN10	Catalyst	10g	
AN20	Catalyst	20g	
AN30	Catalyst	30g	
AN400	Catalyst	400g	

Stainless Steel Anaerobic Jar ANJ Series

ANJ210-7 ANJ240-10 ANJ270-15



Technical characteristics of stainless steel anaerobic culture tank:



◆ Material: stainless steel(SUS304), sturdy and durable, good thermal conductivity;



◆ Temperature resistance: it can withstand high temperature sterilization at 121°C;



◆ Seal: stainless steel lock seal, good sealing performance;



◆ Number of lock catch: 3 pieces.

Specifications

Model	Capacity	Material	Applicable Specification	Standard Configuration	Optional accessories
ANJ210-7	Φ410/h210/7L	SUS304	Place 2 stacks of 12 dishes(Φ9-10cm) or four 250ml conical bottles	Tank body + 2 stainless steel brackets	-
ANJ240-10	Φ240/h240/10L	SUS304	Place 3 stacks of 12 dishes(Φ9-10cm) or five 250ml conical bottles	Tank body + 3 stainless steel brackets	-
ANJ270-15	Φ270/h270/15L	SUS304	Place 4 stacks of 12 dishes(Φ9-10cm) or seven 250ml conical bottles	Tank body + 4 stainless steel brackets	-



Optional accessories

Model	Name	Applicable Specification	Picture
A05077	Connector	1 Suitable for the DWSA05050 connector	
AJ9023	Connector	Suitable for the MARKI connector	
AY6907	Pressure gauge	-	
AN10	Catalyst	10g	
AN20	Catalyst	20g	
AN30	Catalyst	30g	
AN400	Catalyst	400g	