

Protable Leaf Area Meter

BLAM-1 BLAM-2



Technical parameters

Test unit: millimeter, square centimeterPrecision: ±2%
Resolution: 0.01cm2
Measuring Length: ≤1000mm
Measuring Width: ≤160mm
Thickness: ≤8mm
Data capacity: 2000 Groups(Leaf area, average area, length, width)
Package Dimension (W*D*H): 410*220*160mm
G.W.: 1.8kg

Characteristics

- Leaf area, average area, length, width and ratio of length and width.
- PC connected.

Model

Model	Function
BLAM-1	Without computer interface
BLAM-2	USB connects meter to PC, data can be converted to excel format can be printed.

Related Products

- Plant Nutrient Analyzer/Chlorophyll Meter/Plant Photosynthesis Meter/Plant respiration Meter/Plant-Culm-strength Meter