



Facebook



WhatsApp



Twitter



LinkedIn



Instagram



YouTube



INFITEK

LAB & MED equipment solution provider

Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: www.infitek.com

Email: info@infitek.com

Service: support@infitek.com

Address: Room 201, Building A, No.1 Qianwan 1st Road, Qianhai
Shenzhen-Hong Kong Cooperation Zone, Shenzhen (Settled in
Shenzhen Qianhai Business Secretary Co., Ltd.)

Infitek Inc.

522w Riverside Ave Ste N,Spokane,wa 99201

Email: Info@infitek.com



Website

WeChat

INFITEK CO., LTD.

PRODUCT BROCHURE

Slide Cabinet

SLDC-102

Features

- The base is made of SLDC-102 cold-roll steel sheets, the main body and drawers are made of SLDC- 0.8 cold-roll steel sheets.
- This cabinet includes a base and three modules on the top that contain 18 drawers (6 drawers in each module) that accommodate up to 11,000 tissue blocks.
- Movable dividing boards inside the drawers make filing more flexible.
- The hidden handles are made of chrome-coated ABS.
- The slots are specifically manufactured using a metal mold.
- The labeling slots are manufactured through an integrated stamping and forming operation.



Dimensions:

450*478*1295mm (W*D*H)

Slide Cabinet

SLDC-103

Features

- The base is made of SLDC 1.2 cold-roll steel sheets, the main body and drawers are made of SLDC 0.8 cold-roll steel sheets.
- This cabinet includes a base and three modules on the top that contains 18 drawers (6 drawers in each module) that accommodate up to 8928 microscope slides.
- Numbered labels inside the drawers make filing more convenient.
- The hidden handles are made of ABS.
- The slots are specifically manufactured using a metal mold.
- The labeling slots are manufactured through an integrated stamping and forming operation.



Dimensions:

450*478*1295mm (W*D*H)

Slide Cabinet

SLDC-104



Dimensions:
450*315*1210mm (W×D×H)

Features



The base is made of SLDC 1.2 cold-roll steel sheets, the main body and drawers are made of SLDC 0.8 cold-roll steel sheets.



This cabinet includes a base and four stacked modules, a total of 96 slide drying plates can be accommodated, a single module can accommodate 24 slide drying plates, and each compartment can accommodate 12 slide drying plates. The numbered labels in each section make filing more convenient.



The hidden handles are made of ABS.



The slots are specifically manufactured using a metal mold.



The labeling slots are manufactured through an integrated stamping and forming operation.