

# Tissue Water Bath

TWB-1125




## Specification

|                              |                                                                  |
|------------------------------|------------------------------------------------------------------|
| Model                        | TWB-1125                                                         |
| Temperature Adjustment Range | 30~90 °C continuously adjustable, automatic constant temperature |
| Working Temperature          | 30~70 °C                                                         |
| Temperature Control Accuracy | ±2 °C                                                            |
| Ambient Temperature          | 0~40 °C                                                          |
| Electricity                  | AC 220V, 50/60Hz; AC 110V, 50/60Hz                               |
| Power                        | 250W                                                             |
| External Dimension           | 340*340*120mm                                                    |
| Dimension of Cylinder        | 265*190*60mm                                                     |
| Net Weight                   | 4.7kg                                                            |

## Features



Adopts new type heating element and heating block, which has long life and energy saving.



The measured temperature and set temperature are displayed separately.



Made of special black material, corrosion-resistant and wear-resistant.



It has a memory function and automatically retains the set temperature after operation.

# Tissue Water Bath

TSB-1



## Advantages

- High thermal conductivity separate pot design, can quickly heat up and easy to take out water change and cleaning;
- There is an independent high precision temperature controller in the spreading pan which can directly detect the water temperature and accurately and stably control the spreading temperature;
- The water divider can be configured in the pan to support the flexible partition of the water in the pan, and can effectively clean the water debris.



## Features

- Touch interface.
- Intelligent PID temperature control system.
- Triple temperature control protection system.
- Universal angle lighting system.
- Free division function of water area.
- Surface debris cleaning function.
- Separate structure design of water bath dish.

## Specification

| Model                         | TSB-1                       |
|-------------------------------|-----------------------------|
| Dimensions of Water Bath Dish | 247*182*55 mm               |
| Capacity of Water Bath Dish   | 2.3L                        |
| Temperature Range             | 0~99 °C                     |
| Temperature Control Precision | ±1 °C                       |
| Lighting Brightness           | 0~5 level (adjustable)      |
| Working Voltage               | AC 110V~120V / AC 220V~240V |
| Frequency                     | 50Hz / 60Hz                 |
| Power                         | 900W                        |
| Fuse                          | F6A L 250V / F12A L 125V    |
| Safety of Electric            | I                           |
| Ambient Temperature           | 5 °C ~45 °C                 |
| Relative Humidity             | ≤80%                        |
| Air Pressure Range            | 700hPa~1060hPa              |
| External Dimension            | 330*375*130mm               |
| Net Weight                    | 5kg                         |