

Infitek

Hematology Analyzer for Vet

PART AUTO & DIFF FULLY AUTO



Infitek



Infitek



Infitek



Infitek



Infitek

Infitek

Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: infitek.com

Email: info@infitek.com

Service: support@infitek.com

Address: Room 201, Building A, No.1 Qianwan 1st Road, Qianhai
Shenzhen-Hong Kong Cooperation Zone, Shenzhen (Settled in
Shenzhen Qianhai Business Secretary Co., Ltd.)

US Office

INFITEK INC.


522W RIVERSIDE AVE STE N,SPOKANE,WA 99201

EMAIL: INFO@INFITEK.COM



Veterinary Hematology Analyzer, 3Parts

HEMA-V6190



Throughput:60 tests/hour



Sample Volume:10 μ L



21 detection parameters + 3 histograms



Blood Sample:Whole blood, peripheral blood

Features

Rich product performance

- 3-Parts Veterinary Hematology Analyzer, 21 detection parameters + 3 histograms
- Only 10 μ L of venous blood and peripheral blood
- Throughput:60 tests/hour
- Intelligent WBC, RBC, PLT floating landmarks
- Automatic storage, including sample test results, graphics, sample information, and sample abnormality prompts.

Pleasant usage experience

- Multi-user permission management
- Thermal printer, support external printer
- A variety of printing formats are optional, and support customized printing templates
- Various reference value ranges, support user-defined
- Rich quality control calibration function
- English prompt sample abnormal information
- 24-hour power on and customized sleep time

Description

- Veterinary Hematology Analyzer has rich product performance, convenient human interaction and powerful data management system, which can detect different animal samples such as cats, dogs, rabbits, etc.

Convenient human-computer interaction

- 10.4" color touch LCD display, all results and histograms are displayed on the same screen
- English operation, support multiple languages
- Keyboard and mouse input, support barcode scanning (optional)
- One-key maintenance to eliminate faults

Intelligent reagent management function

- Reagents for animals
- Original supporting reagents and maintenance reagents
- Graphic and text display of reagent information and remaining amount
- Supporting consumables: high, medium and low quality control blood, concentrated cleaning solution

Powerful data management system

- Wireless connection: WiFi, Bluetooth and other wireless data transmission, Enables wireless printing, LIS connectivity, lab LAN
- Abundant interfaces: USB, LAN, RS232, HDMI
- Bidirectional LIS communication, HL7 and other transmission formats
- Support remote diagnosis, data import and export, offline reading

After-sales service

- The online after-sales service order system realizes the online contactless order to change the traditional order mode.
- Accurate and convenient after-sale installation, repair and maintenance services.

Specification

Model	HEMA-V6190
21 Parameters	WBC, LYM#, MID#, GR#, LYM%, MID%, GR%, RBC,HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT, P-LCC, P-LCR
3 Histograms	WBC,RBC, PLT
Species	12 pre-set species including Canine, Feline, Equine, Bovine, Murine, Swine, Leporid 4 optional species
Blood Sample	Whole blood, peripheral blood
Sample Volume	10µL
Storage Capacity	Auto-save, more than 200,000 analysis results available
Throughout	60 tests/hour
Data Search Mode	By Sample ID, Patient ID, Date, Multi-keywords. Loading/extracting data, offline reading
Display	10.4 inch colorful touch screen
Operation Mode	Touch screen, keyboard/mouse
Quality Control	L-J, X-R, X-B,etc
Printout	Thermal printer (defaulted),available for external printers Optional printing templates and format
Data Communication	LIS, HL7
Wireless Connection	Wifi, Bluetooth
Ports	LAN×1, RS232×1, HDMI×1
Optional Accessories	Handheld Barcode Scanner
Electricity	110-240V, 50/60Hz
N.W./G.W.(kg)	17/31
Dimension(W*D*H)	440*350*440mm
Package Dimension (W*D*H)(mm)	680*590*650