

Infitek

WATER DISTILLER



Infitek



Infitek



Infitek



Infitek



Infitek

Infitek

Infitek Co., Ltd.

TEL: +86-531-88982330

FAX: +86-531-88983691

Website: infitek.com

Email: info@infitek.com

Service: support@infitek.com

Address: Room 201, Building A, No.1 Qianwan 1st Road, Qianhai
Shenzhen-Hong Kong Cooperation Zone, Shenzhen (Settled in
Shenzhen Qianhai Business Secretary Co., Ltd.)

US Office

INFITEK INC.

522W RIVERSIDE AVE STE N,SPOKANE,WA 99201

EMAIL: INFO@INFITEK.COM



Water Distiller

WDST-S5 WDST-S10

Features

- Easy operation, stable water output.
- Energy saving and emission reduction effects.
- Anti-scalding shell to avoid burns during operation.
- There are two installation modes: desktop or wall mounted.
- Made of 304 high-quality stainless steel, light and beautiful in appearance.
- Equipped with anti-dry burning protection device to prevent dry burning and heating up without water.



Description

The whole machine is made of high quality stainless steel. The separate design structure of the controller and the body reflects features such as easy operation, heat insulation, and safe application. Meanwhile, the built-in distillation shunt structure ensures the stable water output.

Specifications

Model	WDST-S5	WDST-S10
Configuration	Semi-auto	Semi-auto
Material	Stainless steel (304)	Stainless steel (304)
Water Output	5L/h	10L/h
Water Category	Regular water	Regular water
Conductivity (25°C)	<0.5mS/m	<0.5mS/m
Electricity	220V/3kW	380V/7.2kW
External Dimension	310*250*320mm	370*270*370mm
N.W./G.W.	7.5kg/9kg	10.5kg/12kg
Package Dimension	530*340*400mm	580*340*420mm