

# Paraffin Dispenser

PRDP-BMR



## Features



Utilizes full computer automated control and intelligent digital temperature control;  
LED displays the heating status



Utilizes new type of heating elements that enable fast,energy efficient and reliable heating



Temperature sensors are made in USA, precise and reliable;



Memory and automatic restoration function enables the preset temperature to be automatically saved after device operation



Automatically saves the latest operation settings

## Specifications

Model	PRDP-BMR
Temperature Range	Continuously adjustable within 0 - 100 ℃
Temperature control precision	±1℃
Capacity	10000 ml
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	1200W
Dimensions	355×410×540mm (W×D×H)
Net weight	17kg
Package Dimension (W*D*H) (mm)	500*470*570
G.W.(kg)	20kg