

Flocculation Tester

FT-M4A FT-M6A FT-M4C FT-M6C

- Speed: 10~450rpm
- Timer: 0~999 minutes
- Utilizes permanent magnet DC motors
- Equipped with digital speed display and timing functions

Application

It is widely used in universities, environmental protection agencies, CDCs, and the pharmaceutical, metallurgical, chemical, and food industries.

Features

- Utilizes permanent magnet DC motors to ensure high torque, stable rotation, and quiet operation. The speed is fully adjustable, and the system is designed for robust, continuous use.
- The stirring paddles are made of stainless steel, providing high resistance to corrosion.
- This lifting platform is adjustable and very convenient to operate.
- Equipped with digital speed display and timing functions.



Specification

Model	FT-M4A	FT-M6A	FT-M4C	FT-M6C
Type	Asynchronous	Asynchronous	Synchronous	Synchronous
Stirring Places	4	6	4	6
Speed	10~450rpm	10~450rpm	10~450rpm	10~450rpm
Timer	0~999 minutes	0~999 minutes	0~999 minutes	0~999 minutes
Electricity	220V/50Hz	220V/50Hz	220V/50Hz	220V/50Hz
Power of Motor	100W*4	100W*6	260W	260W
Motor	Four motors drive four stirring heads	Six motors drive six stirring heads	One motor drives four stirring heads	One motor drives six stirring heads